

Imagine... consistently superior
cytology preparation



 **Cyto-Tek® 2500**
SAKURA
Cytocentrifuge

CYTO-TEK® 2500 CYTOCENTRIFUGE

The Cyto-Tek® 2500 Centrifuge provides optimal cell recovery and advanced preservation of cellular structure in non-gynecological monolayer slide preparation.

The 2500 was designed with enhanced features focusing on user safety, including the introduction of a new sealable rotor housing with lid which significantly reduces aerosol exposure.

The housing is completely removable and transportable, providing added convenience.

Technicians now have the ability to load and unload specimens in vented areas away from the centrifuge, under a fume hood.

The Cyto-Tek® 2500 Centrifuge has a sleek and modern design and provides the superior clinical performance users have come to expect from Sakura's Cyto-Tek® Centrifuges.

Dimensions

W x D x H 375 x 490 x 225 mm

Weight 11.7 kg


Features

Capacity 12

RPM 200 ~ 2.500 by 10rpm increments

Program Memory 30

Autoclavable Yes (Guard Bowl)



ITEMS & OPTIONS

4300	Cyto-Tek®2500, 100-240V; 50/60 Hz
4301	Chamber Kit 1mL; 50/pcs
4326	Cyto-Tek®6/12 mL Base Holder; 12/pcs
4327	Cyto-Tek®6 mL Gasket; 25/pcs
4328	Cyto-Tek®12 mL Fluid Chamber; 12/pcs
4329	Cyto-Tek®1 mL Fluid Chamber; 200/pcs
4331	Cyto-Tek®6 mL Fluid Chamber; 12/pcs
4333	Cyto-Tek®1 mL Filter Paper; 200/pcs
4334	Cyto-Tek®6/12 mL Chamber Cap; 25/pcs
4335	Cyto-Tek®1 mL Chamber Cap; 200/pcs
4336	Cyto-Tek®1 mL Base Holder; 200/pcs
4337	Cyto-Tek®12 mL Gasket; 25/pcs

F60-900-00	Secondary box for Cyto-Tek®2500, 230V
4338	Cyto-Tek®6 mL Filter Paper; 500/pcs
4339	Cyto-Tek®12 mL Filter Paper; 500/pcs
8704	Giemsa; 2.5 L
8705	Giemsa; 5 L
8716	Papanicolaou 2A OG6; 2.5 L
8717	Papanicolaou 2A OG6; 5 L
8719	Papanicolaou 2B Orange II; 2.5 L
8720	Papanicolaou 2B Orange II; 5 L
8722	Papanicolaou 3B EA 50; 2.5 L
8723	Papanicolaou 3B EA 50; 5 L
8731	Cervix Spray Fixative; 125 mL



First we understand.
Then we innovate.



Sakura Finetek Europe B.V.
The Netherlands
Tel: +31 (0)88 592 00 00
Fax: +31 (0)88 592 00 01
Sakura@sakura.eu
www.sakura.eu